

# CELL CULTURE CONTAMINATION

InvivoGen *infocus*

PRACTICAL GUIDE



## OVERVIEW OF MICROBIAL CONTAMINANTS

Whether remaining unnoticed or expanding rapidly, microbes can seriously alter cell morphology and functions, becoming a serious threat to your research. InvivoGen offers microbial detection tools and elimination reagents, as well as preventive tips.

- Protect your cells
- Mycoplasma contaminations
- Bacterial contaminations
- Endotoxin contaminations
- Fungal contaminations

# InvivoGen *infocus*

## A PRACTICAL GUIDE TO AVOID CELL CULTURE CONTAMINATION

**4**  
PAGE

### ***YOUR CELLS ARE PRECIOUS, PROTECT THEM!***

Detection  
Prevention  
Elimination

**5**  
PAGE

### ***MYCOPLASMA CONTAMINATIONS***

Mycoplasma features & dangers  
Detection of mycoplasma contamination  
Elimination of mycoplasma

**10**  
PAGE

### ***BACTERIAL CONTAMINATIONS***

The usual suspects  
How to detect bacterial contamination?  
Elimination of bacteria

**12**  
PAGE

### ***ENDOTOXIN CONTAMINATIONS***

What are endotoxin features?  
What are the risks for my experiments?  
Detection of endotoxins in biological reagents

**14**  
PAGE

### ***FUNGAL CONTAMINATIONS***

How to detect fungal contamination in cell cultures?  
Elimination of fungi

**15**  
PAGE

### ***SUMMARY TABLE OF INVIVOGEN'S ANTIMICROBIAL AGENTS***

### ***REFERENCES***


**M**icrobial contamination of cell cultures is a serious and relentless threat to your research. Invasive mycoplasma, bacteria, and fungi can kill or drastically alter cells in culture, leading to disastrous results, lost time, and wasted resources. This brochure provides an insight into the contaminants that are most likely to invade your cultures, the good practices to avoid them, and the solutions to eliminate them. As experts in innate immunity and microbiology, we know how these biological contaminants can interfere with experimental results. While bacterial and fungal contaminations are eventually detected by the naked eye, mycoplasma and endotoxins remain invisible. Undetected contaminants are a serious concern, as they may have led to data misinterpretation, many of which have been published. As a consequence, journals now frequently ask for evidence of absence of mycoplasma and endotoxins in cell cultures. Moreover, pharmaceutical companies developing future therapeutics cannot afford contaminations as they will compromise their research and reputation. At InvivoGen, we strive for excellence. We provide high-quality products and mycoplasma-free cell lines all around the world. This guide will help you address every stage of microbial infection, and choose the right InvivoGen product to detect, eliminate, and prevent contaminations in your cell cultures.

## ABBREVIATIONS

**CFU**: colony-forming unit  
**DNA**: deoxyribonucleic acid  
**IRF**: interferon regulatory factor  
**kDa**: kilodalton  
**LAL**: limulus amoebocyte lysate  
**LPS**: lipopolysaccharide

**NF-κB**: nuclear factor kappa-light-chain enhancer of activated B cells  
**PCR**: polymerase chain reaction  
**PRR**: pattern recognition receptor  
**rRNA**: ribosomal ribonucleic acid  
**SEAP**: secreted embryonic alkaline phosphatase  
**TLRs**: toll-like receptors

# TINY ORGANISMS, BIG HEADACHES



1956

First detection and isolation of mycoplasma in cell cultures

+180

Different mycoplasma species have been identified

AT LEAST 20


Distinct mycoplasma species isolated from contaminated cell cultures

95%

of cell culture contamination is due to 6 mycoplasma species

- *M. orale*
- *M. fermentans*
- *M. hyorhinis*
- *M. hominis*
- *M. arginini*
- *A. laidlawii*

## MYCOPLASMA IN NUMBERS



80% OF LAB OPERATORS CARRY MYCOPLASMA



## CONTAMINATION PREVENTION TIPS

Quarantine all new cell cultures and animal products entering the laboratory

Be vigilant and always practice good aseptic technique

Avoid talking over your cells  
6% of operators spread mycoplasma by talking

Avoid sneezing  
38% of operators spread mycoplasma by sneezing

Routinely test your cells

Ask manufacturers for certification proving the absence of mycoplasma contamination

InvivoGen

OFFERS

>14

ANTI-MICROBIAL SPECIFIC REAGENTS & ASSAYS TO ACCOMPANY YOUR SUCCESSFUL RESEARCH

Mycoplasma contamination testing status must also be reported.

section. Editors reserve the right to demand that the data be removed from the paper if the justification is deemed unsatisfactory. In addition, authors must identify the source of cell lines (with catalog number). Cell lines have been authorized for use in research only. No cell line authentication testing will be provided upon request. As of May 2015, there is no issue of cell line authentication. It is good practice to obtain cell lines from reputable repositories, to routinely authenticate cell line stocks and test them for mycoplasma contamination. Resources on cell line authentication follow.

It is good to obtain cell lines from reputable repositories, to routinely authenticate cell line stocks and test them for mycoplasma contamination

# YOUR CELLS ARE PRECIOUS, PROTECT THEM!

## Detection

Microbial contamination must be detected **as early as possible**. Detection methods depend on the nature of the microbe. They include **biological assays**, **PCR**, **fluorescence** or **chemical staining**, **optical microscopy**, **turbidimetry**, **pH** measurements, or simple **visual inspection**. Bacteria and fungi can usually be identified by optical microscopy. Their **fast growth rate** allows their detection by the naked eye as early as 48 hours (i.e. over the weekend), the contaminated cultures appearing turbid or spotty. Subsequently, identification of these micro-organisms can be performed with testing kits. Mycoplasma in cell cultures cannot be detected visually, not even by optical microscopy. Hence, these microbes can go unnoticed for long periods and are identified using dedicated assays.

## Prevention

Knowing the sources of microbial contamination is crucial for **minimizing the risk to cell cultures** (see below). Although absolute prevention is impossible, you can take various measures to prevent infection. Firstly, ensure that you are working in a **sterile environment** and using proper **aseptic technique**. Secondly, **quarantine** any incoming cell cultures until these have been confirmed free of contamination. Thirdly, **monitor** your cell cultures for contamination on a regular basis by optical microscopy and detection kits. Lastly, you can use **antibiotic cocktails**, such as those offered by InvivoGen, specifically designed for taking a preventive strike against microbes that would be difficult to detect in new cultures (i.e. primary cells, or cloning).

## Elimination

Typically, once invasive microbes are detected in cell cultures it is recommended to discard the cells and the media. However, some cell cultures are so precious that they cannot be lost (i.e. stable clone selection, cell lines derived from explanted tissues, primary cells) and are not available elsewhere (i.e. not yet frozen). In such situations, InvivoGen provides antibiotics to eradicate the contamination **surely** and **rapidly without damaging** your cells.

### 1. INFRASTRUCTURE:

Fume-hoods, ventilation units and laboratory furniture can house surface microorganisms and spread airborne ones. Clean working areas daily with alcohol and monthly with bleach. Regularly change all air filters, and empty all media traps at least weekly.

### 2. OPERATOR:

Laboratory staff can transmit microbes and dust from their skin, clothes and bodies to cell cultures. Wear proper safety garments and use aseptic technique.

### 3. PIPETTES, TIPS, SYRINGES & VACUUM PUMPS:

A contaminated pipette can destroy multiple cell cultures. Use sterile pipettes, tips and syringes and never reuse disposables. Make sure to empty and clean the vacuum pump reservoir and tubings regularly.

### 4. EQUIPMENT:

Glassware, incubators and water baths can easily be contaminated. Keep all equipment sterile and frequently change the water in the baths, which are notorious for fungal infections. Regularly disinfect all incubators.

### 5. SAFETY GARMENTS:

Labcoats and gloves must be clean and should be worn only in the lab. Use disposable garments only once and discard them immediately afterwards.

### 6. CULTURE MEDIA, REAGENTS & CELLS:

Before using any culture medium, reagent or cells, confirm that they are sterile. Carefully seal all containers.






# MYCOPLASMA CONTAMINATIONS

**M**ycoplasma are the smallest and simplest self-replicating organisms. Because of their small size (100 nm) and lack of a rigid cell wall, mycoplasma are undetectable by visual inspection, pass through standard filtration, and are resistant to a large number of antibiotics<sup>1</sup>. Mycoplasma contamination is a major problem in cell culture, affecting the validity of experimental results as well as the quality and safety of cell-based biopharmaceuticals<sup>2</sup>.

## Mycoplasma features



**Mycoplasma** belong to the class of Mollicutes, which members are distinguished by their lack of cell wall and their plasma-like form. Mycoplasma are **highly infectious**, for all types of eukaryotic cells, including primary cells. Hundreds of mycoplasma can attach to a single cell, fuse with the cell membrane, multiply, and eventually outnumber cultured cells by 1000-fold. Mycoplasma can drastically **alter cell cultures** and **skew research results** (see below). Mycoplasma lipoproteins are potent activators of immune cells upon sensing by **Toll-like receptor 2** (TLR2), their preferential pattern recognition receptor (PRR)<sup>3</sup>. The absence of a cell wall in mycoplasma confers them resistance to commonly used antibiotics, such as penicillin and streptomycin. Moreover, their **tiny size (~100 nm)** does not allow their elimination by standard 0.2 µm filtration. Thus, major precautions are required to prevent contamination of cell cultures.


All InvivoGen's cell lines are guaranteed mycoplasma-free based on Plasmotest™ results

## Major impacts of mycoplasma contamination on cell functions

Mycoplasma **compete** with host cells for nutrients and biochemical precursors. As a consequence, they alter many cell functions, such as **cell metabolism** and **cell growth**, ultimately leading to **cell death**. A microarray analysis on contaminated cultured human cells has revealed the severe effects that mycoplasma can have on the expression of hundreds of genes, including some that encode receptors, ion channels, growth factors, and oncogenes<sup>4</sup>. Upon adhesion or fusion with the host cell membrane, they can cause further damage to the cell by interfering with **signaling cascades** and **cytokine production**<sup>5</sup>. These detrimental effects can strongly impact scientific results and invalidate the findings of a study, especially when it involves immune cells expressing TLR2, such as macrophages<sup>3,6</sup>.

### TOP 5 REASONS TO TEST

- ➔ Loss of precious cell lines
- ➔ Serious impact on data reliability and reproducibility
- ➔ Testing is required by most journals for publication
- ➔ Up to 35% worldwide contamination, including cell banks
- ➔ Loss of time and money



# DETECTION OF MYCOPLASMA CONTAMINATION



InvivoGen offers two mycoplasma detection kits to allow timely intervention. Both kits allow fast and accurate detection of mycoplasma species that most commonly contaminate cell culture. MycoStrip™ is based on 'immediate' (~1 h) detection of mycoplasma genomic content using immunochromatographic strips. Plasmotest™ is a colorimetric cellular assay based on detection of mycoplasma lipoproteins.

Catch mycoplasma off-guard using 'immediate' and 'routine' tests

## Genomic detection strips

### MycoStrip™

NEW

- **Recommended for immediate results**
- **Simple:** No special lab equipment required
- **Rapid:** Hands-on time <15 min. Total duration: 1 h.
- **Clear:** One band – negative for mycoplasma  
Two bands – positive for mycoplasma

Detection of cell culture contaminating mycoplasma by MycoStrip™ is based on **isothermal PCR**. The **16S rRNA gene** for the most commonly found mycoplasma species in cell culture, accounting for 95% of contamination, is targeted and amplified using our proprietary Reaction mix. Results are visualized as a band on an **immunochromatographic strip** within 5 minutes.

PRODUCT	DESCRIPTION	QTY	CAT. CODE
MycoStrip™	Mycoplasma contamination detection kit (strips)	10 tests	rep-mys-10
		20 tests	rep-mys-20
		50 tests	rep-mys-50
		100 tests*	rep-mys-100

\*without cassette

## FREQUENTLY ASKED QUESTIONS



### Q: Do MycoStrip™ and Plasmotest™ detect only live mycoplasma?

A: Both kits detect live and dead mycoplasma. MycoStrip™ and Plasmotest™ detect the presence of mycoplasma DNA and lipoproteins, respectively.


### Q: Does Plasmotest™ detect only mycoplasma in cell cultures?

A: Plasmotest™ relies on the activation of TLR2. Therefore, it can detect both mycoplasma and bacteria contaminants. However, while the mycoplasma cannot be detected by the naked eye, bacteria contamination is visible and leads to a decreased pH (change of medium color to yellow) and medium turbidity due to bacterial growth.

More FAQs online

[www.invivogen/mycostrip](http://www.invivogen/mycostrip) & [www.invivogen/plasmotest](http://www.invivogen/plasmotest)

Mycoplasma & Acholeplasma




Isothermal PCR




Immunochromatographic strip

### MycoStrip™ procedure


01 Combine: Reaction mix + Reaction buffer + Sample



02 Perform isothermal PCR



03 Add Migration buffer and transfer sample onto the cassette



04 Wait 2 to 5 minutes for results



## Colorimetric cellular assay


### PlasmoTest™

- **Recommended for routine tests**
- **Simple:** Colorimetric detection in cell culture supernatant
- **Rapid:** Hands-on time <1 h. Results overnight.
- **Sensitive:** LOD  $5 \cdot 10^2$  -  $5 \cdot 10^5$  CFU/ml of culture supernatant


PlasmoTest™ relies on Toll-like receptor 2 (TLR2), the preferential pattern recognition receptor (PRR) for **mycoplasma lipoproteins**<sup>1</sup>. Our proprietary **HEK-Blue™-2 cells** stably express human TLR2 and an NF-κB-inducible SEAP (secreted embryonic alkaline phosphatase) reporter gene. They are cultured in HEK-Blue™ detection medium, make them ideal for routine tests. **The simple addition of test samples** to these cells provides **colorimetric results** with sensitivity similar to luminescence-based biochemical assays. SEAP activity can be measured at 620-655 nm using a spectrophotometer. The absorbance is in direct proportion to the amount of contaminants. For your convenience, HEK-Blue™-2 cells are provided with Normocin™, an antibiotic cocktail to prevent cell culture contamination with mycoplasma, bacteria, and fungi (see page 11).

Mycoplasma & Acholeplasma

Lipoproteins



Reporter cell line




Purple/blue: Positive

Pink: Negative


Colorimetric read-out

### PlasmoTest™ procedure


#### 01 Inactivate endogenous alkaline phosphatase in the sample



#### 02 Add heated sample to HEK-Blue™-2 sensor cells cultured in HEK-Blue™ detection medium




#### 03 Mycoplasma lipopeptide detection by HEK-Blue™-2 cells



#### 04 SEAP reporter activity

- Naked-eye: Purple/blue: Positive
- Pink: Negative
- Measure OD 620-655nm (optional)



### MYCOPLASMA DETECTION METHODS



	MycoStrip™	PlasmoTest™
Target	16S rRNA gene	Lipoproteins
Method	Isothermal PCR	Reporter bioassay
Ease of use	++++	++
Specificity	At least the 8 most common species	
Sensitivity	$10 \cdot 10^2$ CFU/ml	$10^2 \cdot 10^5$ CFU/ml
Experiment duration	<1 hr	<1 hr (hands on) OVN (incubation)
Additional equipments/ reagents	Heating bath/block	Incubator Spectrophotometer
Naked-eye data visualization	Yes	Yes
Cost	Reasonable	Cost-effective

### WHAT IF MY TEST IS POSITIVE ?

**Your culture is easily treatable with InvivoGen's anti-mycoplasma reagents.**

Treat your culture and eradicate the contamination using **Plasmocin™** or **Plasmocure™**. Upon completion of the treatment (~2 weeks), re-test using MycoStrip™ or PlasmoTest™ comparing your newly treated culture with your previous sample.

PRODUCT	DESCRIPTION	QTY	CAT. CODE
PlasmoTest™	Mycoplasma contamination detection kit (cells)	1 kit (250 tests)	rep-pt1
PlasmoTest™ Controls	Controls for PlasmoTest™ detection kit	200 tests	pt-ctr2
PlasmoTest™ Refills	Reagents for PlasmoTest™ detection kit	500 samples	rep-ptrk

# ELIMINATION OF MYCOPLASMA



InvivoGen has over 40 years of experience in developing anti-mycoplasma solutions for the scientific community. Plasmocin™ and Plasmocure™ are unique anti-mycoplasma reagents that combine two distinct sets of antibiotics in a single ready-to-use product. They act fast with little to no cytotoxicity, ensuring you save precious cell lines and data.

Plasmocin™ and Plasmocure™ allow broad and rapid elimination of mycoplasma contaminants.

## A preventive & removal treatment

### Plasmocin™

- **Recommended for broad mycoplasma elimination**
- **Fast:** Rescues cell cultures in 2 weeks
- **Safe:** Little to no toxicity on mammalian cells
- **2 formats:** Prophylactic and Treatment

Plasmocin™, a frequently cited mycoplasma removal agent<sup>7-11</sup>, is effective against [all common mycoplasma strains](#), both extracellular and intracellular. For maximum efficiency, Plasmocin™ contains a formulation of two antibiotics: the first one blocks [protein synthesis](#), and the second one stops [DNA replication](#). A component of Plasmocin™ is actively transported into mammalian cells, ensuring that following treatment, cell cultures do not become re-infected. Thus, Plasmocin™ is more effective than other reagents on the market in eradicating mycoplasma and preventing resistant strain generation<sup>7</sup>.

- Plasmocin™ **prophylactic** is a product that can be used on a [regular basis](#) to [prevent](#) mycoplasma contaminations.
- Plasmocin™ **treatment** is intended for mycoplasma [elimination](#) within [2 weeks](#).

## An alternative removal treatment

### Plasmocure™

- **Recommended for Plasmocin™-resistant strains**
- **Fast:** Rescues cell cultures in 2 weeks
- **Safe:** Little to no toxicity on mammalian cells
- **Reliable:** Extremely low mycoplasma regrowth rate


Plasmocure™ is a [second-line anti-mycoplasma reagent](#) that potently eradicates Plasmocin™-resistant mycoplasma<sup>12</sup>. It combines two antibiotics that act through different mechanisms than those found in Plasmocin™. The first antibiotic binds to the 50S subunit of the ribosome and blocks [peptidyltransferase activity](#). The second antibiotic binds to isoleucyl-tRNA synthetase and halts [isoleucine incorporation](#) into mycoplasma proteins. A two-week treatment with Plasmocure™ is usually sufficient to completely eliminate mycoplasma. If mycoplasma elimination is not completed after two weeks, Plasmocure™ can be administered for an additional week. A moderate slowdown of cell growth may be observed, but full recovery of the cell line is expected once mycoplasma are eliminated.

## Tips for successful mycoplasma elimination



### OPTIMAL TREATMENT CONCENTRATION

We recommend to test 3 different treatment concentrations according to the protocol, including a no-treatment condition.



### BACK-UP CELLS

Maintain a culture duplicate and/or frozen vial without treatment:

- to validate the treatment efficacy
- to start over with new treatment conditions if necessary.



### ALTERNATIVE TREATMENT

In case of a lack of Plasmocin™ efficacy, use Plasmocure™ (or vice versa).




### REGULAR TESTING

Early mycoplasma detection, using MycoStrip™ or Plasmotest™, allows timely intervention. We recommend testing cell cultures every ~2-3 weeks.



## Mycoplasma surveillance and elimination

Upon reception of a new cell line or maintenance of cell cultures



Use **Plasmocure™** for **Plasmocin™**-resistant mycoplasma contaminations

### **DON'T STRESS**

**Many infected cell lines have been successfully treated with Plasmocin™**

Including cancer cell lines<sup>9</sup>, virus-producing cells<sup>10</sup>, induced pluripotent stem cells<sup>11</sup>, and human embryonic stem cells<sup>13</sup> with no permanent alterations.

## MYCOPLASMA ELIMINATION REAGENTS



	Plasmocin™	Plasmocure™
Treatment duration	2 weeks	2 weeks
Ease of use	++++	++++
Efficacy	+++	++++
Cytotoxicity	+/-	+/-
Resistance <sup>11</sup>	-	--

## FREQUENTLY ASKED QUESTIONS



**Q: Can I use Plasmocin™ during the initial culture phase before making my frozen stocks?**

A: Yes, it is even recommended. However, it is necessary to check that the cells are not contaminated using MycoStrip™ or Plasmotest™.

**Q: What is the toxicity of Plasmocin™ and Plasmocure™ to eukaryotic cells?**

A: Plasmocin™ and Plasmocure™ targets are absent in eukaryotic cells, ensuring low cytotoxicity.

**Q: What to use? Plasmocin™ or Plasmocure™?**

A: If your cells are positive for mycoplasma, we recommend starting a treatment with Plasmocin™. Plasmocure™ should be used in the case of resistance to Plasmocin™.

**More FAQs online**

[www.invivogen/plasmocin](http://www.invivogen/plasmocin) & [www.invivogen/plasmocure](http://www.invivogen/plasmocure)

## They trust InvivoGen

- Kapralov A. A. et al., 2020. *Nat Chem Biol*. Redox lipid reprogramming commands susceptibility of macrophages and microglia to ferroptotic death. DOI: 10.1038/s41589-019-0462-8.
- Kazemiha V. M. et al., 2019. *Cell J*. Effectiveness of Plasmocure™ to eliminate mycoplasma species from contaminated cell cultures: A comparative study versus other antibiotics. DOI: 10.22074/cellj.2019.5996


PRODUCT	DESCRIPTION	QTY	CAT. CODE
Plasmocin™ prophylactic	Reagent for preventing mycoplasma contamination	25 mg (10 x 1 ml)	ant-mpp
Plasmocin™ treatment	Mycoplasma elimination reagent	25 mg (1 x 1 ml) 50 mg (2 x 1 ml)	ant-mpt-1 ant-mpt
Plasmocure™	Mycoplasma elimination reagent	100 mg (1 ml)	ant-pc

# BACTERIAL CONTAMINATIONS

Bacteria are a large and ubiquitous group of unicellular micro-organisms. They are typically a few micrometers in diameter and can have a variety of shapes, ranging from spheres to rods and spirals. Bacteria are the most commonly encountered biological contaminants in cell culture<sup>14, 15</sup>. Despite being detectable using a light microscope, bacteria can easily be mistaken for cellular debris, especially during the early stages of contamination.



## The usual bacterial suspects



Bacteria can easily be mistaken for cellular debris, especially at the early stages of contamination.

## Sources of contamination




### LAB OPERATORS

Major source of contamination with *Staphylococcus* species.



### DIRTY WATER BATHS, INCUBATORS, GLASSWARE, ...

Major source of contamination with *Pseudomonas* and *Flavobacterium* species.



### CELLS ISOLATED FROM ANIMAL TISSUES

Contamination from the commensal flora and/or subclinical infections.

## Detection of bacterial contamination



In the **early stages** of contamination, bacteria can be mistaken for cellular debris as they are much smaller (1-10 µm) than eukaryotic cells (10-100 µm). Therefore, it is **important to check your cell cultures** under a **light microscope** using phase contrast (100x - 400x).

Do you see an abnormal presence of small black dots, rods, spirals, either alone, in chains, or clusters? Are they motile? If so, your culture is probably contaminated.



Because of their **fast growth rate**, bacteria cause a **change in the culture medium in just 48 hours**, making the contamination clearly visible with the **naked eye** from 10<sup>5</sup> CFU/ml. The culture medium appears cloudy, and if it contains phenol red, a rapid color change from red to yellow indicates a decrease in pH, a consequence of bacteria metabolism. The culture **environment is no longer suitable for eukaryotic cells** leading eventually to their **death**.

These cultures should be discarded, or if irreplaceable, treated with InvivoGen's antibiotic cocktails (see next page).

# ELIMINATION OF BACTERIA

## Contamination preventive reagent

### Normocin™

- **Broad-spectrum:** Kills mycoplasma, bacteria, and fungi
- **Safe:** Little to no toxicity on mammalian cells

Normocin™ is an innovative formulation of three antibiotics active against **mycoplasma**, **bacteria (Gram+ and Gram-)**, and **fungi**, including yeasts. It is widely used and cited as a “**routine addition**” to cell culture media to prevent contamination. Normocin™ can be used in combination with penicillin and streptomycin (Pen-Strep) solutions to broaden the anti-bacterial spectrum. Normocin™ provides maximum protection against microbial contamination with **minimum cytotoxicity**.

#### DID YOU KNOW?

All InvivoGen's cell lines are provided with Normocin™, a broad-spectrum anti-microbial agent, to help you keep your cells safe

## Anti-microbial agent for primary cells

### Primocin™

- **Broad-spectrum:** Kills mycoplasma, bacteria, and fungi
- **Safe:** Little to no toxicity on mammalian cells

Primocin™ contains four compounds, with three of these blocking DNA and protein synthesis in **Gram+ bacteria**, **Gram- bacteria**, and **mycoplasma**. The fourth compound eradicates **fungi**, including yeasts, by disrupting ionic exchange through the cell membrane. Primocin™ has been successfully used with many **primary cells**, including mouse- and human-tumor-derived cell lines, embryonic cells, and induced pluripotent stem cells.

## Multidrug-resistant bacteria removal agent

### Normocure™

- **Ready-to-use:** Add directly into medium bottles or flasks
- **Fast:** Rescues cell cultures in 3 passages
- **Safe:** Little to no toxicity on mammalian cells

Normocure™ is the best weapon to save your valuable cell lines from **Gram+ and Gram- bacteria, especially non-fermenting Gram-bacteria** that are resistant to Pen-Strep and Normocin™. Normocure™ is a cocktail of three components belonging to different antibiotic families. After the first passage, >99% of bacterial contaminants are eliminated. The targets of these antibiotics are absent in eukaryotic cells, ensuring Normocure™'s low cytotoxicity.

## FREQUENTLY ASKED QUESTIONS

**Q: Do the anti-microbial agents you have in your catalog interfere with selective antibiotics?**

A: No, they do not interfere with common selective antibiotics such as G418, Blastcidin, Puromycin, Hygromycin B, or Zeocin™.

**Q: We have a bacterial contamination but cannot determine which bacterial strain has contaminated our cultures. What would be the best option to ensure bacteria elimination?**

A: We would highly recommend using Normocure™ which is a broad-spectrum antibacterial agent highly effective against Gram+ and Gram-bacteria. Cell cultures contaminated with bacteria from the environment, such as *Staphylococcus* species and *Achromobacter* species, can be efficiently cured by Normocure™ treatment.

**More FAQ online**

[www.invivogen/normocin](http://www.invivogen/normocin), [www.invivogen/normocure](http://www.invivogen/normocure), [www.invivogen/primocin](http://www.invivogen/primocin)

## They trust InvivoGen

- Magupalli V. G. *et al.*, 2020. *Science*. HDAC6 mediates an aggresome-like mechanism for NLRP3 and pyrin inflammasome activation. DOI: 10.1126/science.aas8995.
- Polyzoos A. A. *et al.*, 2019. *Cell Metab*. Metabolic Reprogramming in Astrocytes Distinguishes Region-Specific Neuronal Susceptibility in Huntington Mice. DOI: 10.1016/j.cmet.2019.03.004.
- Ganesh K. *et al.*, 2020. *Nat. Cancer*. L1CAM defines the regenerative origin of metastasis-initiating cells in colorectal cancer. DOI: 10.1038/s43018-019-0006-x.

PRODUCT	DESCRIPTION	QTY	CAT. CODE
Normocin™	Reagent for preventing microbial contamination	500 mg (10 x 1 ml) 1 g (1 x 20 ml)	ant-nr-1 ant-nr-2
Normocure™	Microbial contamination elimination reagent	100 mg (2 x 1 ml)	ant-noc
Primocin™	Reagent for preventing microbial contamination in primary cells	500 mg (10 x 1 ml) 1 g (1 x 20 ml)	ant-pm-1 ant-pm-2


# ENDOTOXIN CONTAMINATIONS




**E**ndotoxins, also known as lipopolysaccharides (LPS) or lipoglycans, are a major cell wall component of Gram-negative bacteria. Sources of endotoxins include media, sera, water, buffers, and other cell culture reagents, such as trypsin. Endotoxins are potent inducers of inflammatory responses both *in vitro* and *in vivo*. Extra care needs to be taken with solutions and reagents that are sterile but may still contain bacterial components, and monitoring for the presence of endotoxins in cell culture reagents is crucial.

Endotoxin contaminations are a major concern for cell cultures and production of injectable drugs.

## Endotoxin features




## Sources of contamination



## Risks for *in vitro* and *in vivo* experiments

The lipid A moiety of LPS activates the **Toll-like receptor 4 (TLR4)** at the cell surface or in endosomes, and subsequently induces the activation of **MAP Kinases**, **NF-κB**, and **IRF (interferon regulatory factor)** pathways<sup>20, 21</sup>. Although the effects of endotoxins vary according to concentration and cell type, it has been shown that these molecules can alter **cellular morphology**, **proliferation**, and **transfection efficiency**<sup>22</sup>. Moreover, LPS induces the production of **inflammatory cytokines** such as TNF-α, IL-1β, IL-6, IL-8, and IL-18, and activation of the caspase 4/5/11-NLRP3 non-canonical inflammasome and cause **pyroptotic cell death**<sup>20, 23</sup>.

*In vivo*, low endotoxin concentrations can stimulate the **immune system**, while high concentrations can induce **fever**, **hypotension**, **multiple organ failure**, and even **death**<sup>22</sup>.



# DETECTION OF ENDOTOXINS IN BIOLOGICAL REAGENTS

## Colorimetric cellular assay

### HEK-Blue™ LPS Detection Kit 2

- **Versatile:** Detection in virtually all biological reagents
- **Highly sensitive:** Detects as little as 0.01 EU/ml
- **Economical:** 1 kit allows to perform up to 500 tests


This kit relies on InvivoGen's proprietary **HEK-Blue™-4 cells**, which stably express human **TLR4** and an **NF-κB-inducible secreted embryonic alkaline phosphatase (SEAP) reporter gene**. The simple experimental procedure makes this kit ideal for routine tests. The presence of minute quantities of LPS in test samples incubated overnight with the HEK-Blue™-4 cells leads to the activation of NF-κB and expression of the SEAP reporter. The activity of SEAP in the culture supernatant is assessed by QUANTI-Blue™ Solution, a colorimetric SEAP detection medium, and measured at 620-655 nm using a spectrophotometer. The absorbance is in direct proportion to the amount of endotoxin present. The endotoxin concentration can be calculated from a standard curve obtained using serial dilutions of the HEK-Blue™ Endotoxin Standard provided in the kit.

#### TOP 5 REASONS TO USE



- ➔ **Sustainable alternative to the costly and laborious limulus amoebocyte lysate (LAL) assay.** Freeze HEK-Blue™-4 cells and re-order the other components separately.
- ➔ Unlike the LAL assay, HEK-Blue™-4 cell-based LPS detection **does not lead to false positives** if the sample contains (1,3)-β-D-glucan.
- ➔ **For all biological samples**, including medium, serum, chemical preparations, and vaccine adjuvants.
- ➔ **Convenient visualization of results by change of medium color** from pink to purple/blue.
- ➔ **Includes a booklet for experimental procedures, a graphical method to calculate the endotoxin concentration, and a troubleshooting section.**

#### Principle of LPS Detection Kit 2



#### They trust InvivoGen

- Gray M. A. *et al.*, 2020. *Nat Chem Biol*. Targeted glycan degradation potentiates the anticancer immune response in vivo. DOI: 10.1038/s41589-020-0622-x.
- Imbert P. R. C. *et al.*, 2021. *Curr Biol*. An Acquired and Endogenous Glycocalyx Forms a Bidirectional "Don't Eat" and "Don't Eat Me" Barrier to Phagocytosis. DOI: 10.1016/j.cub.2020.09.082.

PRODUCT	DESCRIPTION	QTY	CAT. CODE
HEK-Blue™ LPS Detection Kit 2	Assay for the detection and quantification of biologically active LPS	1 kit	rep-lps2
HEK-Blue™ Selection	Antibiotics for maintenance of HEK-Blue™ cells	10 x 1 ml	hb-sel
HEK-Blue™ Endotoxin Standard	Standardized <i>E. coli</i> 055:B5 LPS	10 x 50 EU	rep-hbes-10
QUANTI-Blue™ Solution	Alkaline phosphatase detection medium (liquid form, 100X)	5 ml 10 ml	rep-qbs rep-qbs2



# FUNGAL CONTAMINATIONS



Unlike mycoplasma and bacteria, fungi are eukaryotes and can exist as round or oval bodies (yeasts) that can form chains or clusters, or as long thin filaments (hyphae)<sup>24</sup>. Molds are a group of hyphae that appear as fuzzy patches in the advanced stages of contamination<sup>24</sup>. Fungal contamination is a major problem in cell culture, affecting the validity of experimental results. More importantly, these types of contaminations are extremely hard to eradicate, as fungi can spread via spore mobility in the air. Dormant spores of many fungal species can survive in extremely harsh and inhospitable environments, only to become activated when they encounter suitable growth conditions.


Cell cultures can often be cured of fungal contamination when detected early and treated with Fungin™ antimycotic agent.

## How to detect fungal contamination in cell cultures?

**Yeasts** are the smallest form (3-10 µm) of fungi. They can be seen using a light microscope.

**Hyphae** can be detected by the naked eye or a light microscope depending on their size and growth stage.

In the case of substantial contamination, colonies form as **molds** floating on the surface. In this case, do not open the vessel and discard the cultures to avoid spreading spores. Sometimes, the medium pH may increase, resulting in phenol-red-containing media appearing pink.



### They trust InvivoGen

- Busslinger G.A. et al., 2021. *Cell Rep.* Human gastrointestinal epithelia of the esophagus, stomach, and duodenum resolved at single-cell resolution. DOI: 10.1016/j.celrep.2021.108819
- Lebreton F. et al., 2019. *Nat Commun.* Insulin-producing organoids engineered from islet and amniotic epithelial cells to treat diabetes. DOI: 10.1038/s41467-019-12472-3.

## ELIMINATION OF FUNGI

### For preventive and removal treatments

#### Fungin™

- **Effective:** Kills fungi (yeasts, hyphae, and molds)
- **Fast:** Rescues cell cultures in 5-10 days
- **Safe:** Minimum toxicity on mammalian cells
- **Reliable:** Can be used as a prophylactic treatment

Fungin™ is a soluble antimycotic compound that kills the **different forms of fungi (yeasts, hyphae, and molds)** by disrupting ionic exchange through the cell membrane. It is a highly stable compound and does not need to be dissolved in deoxycholate (which is cytotoxic). Therefore, it is an excellent alternative to the use of Amphotericin B antimycotic. Fungin™ can be used at either low or high concentrations for **routine prevention or contamination removal**, respectively. It provides maximum protection against fungal contaminants commonly found in cell culture, such as *Candida albicans* and *Aspergillus spp.* Fungin™ may be added to media containing commonly used antibacterial agents, such as penicillin and streptomycin (Pen-Strep). Fungin™ is a highly referenced antimycotic reagent.

PRODUCT	DESCRIPTION	QTY	CAT. CODE
Fungin™	Reagent for preventing fungal contamination	75 mg (5 x 1.5 ml) 200 mg (1 x 20 ml)	ant-fn-1 ant-fn-2

# SUMMARY TABLE OF INVIVOGEN'S ANTIMICROBIAL AGENTS

CONTAMINANT	PRODUCT	DESCRIPTION	PAGE
Mycoplasma	MycoStrip™	Mycoplasma contamination detection kit	6
	PlasmoTest™	Mycoplasma contamination detection kit	7
	PlasmoTest™ Controls	Controls for PlasmoTest™ detection kit	7
	PlasmoTest™ Refills	Reagents for PlasmoTest™ detection kit	7
	Plasmocin™ prophylactic	Reagent for preventing mycoplasma contamination	8
	Plasmocin™ treatment	Mycoplasma elimination reagent	8
	Plasmocure™	Mycoplasma elimination reagent	8
Bacteria	Normocin™	Reagent for preventing microbial contamination	11
	Normocure™	Microbial contamination elimination reagent	11
	Primocin™	Reagent for preventing microbial contamination in primary cells	11
Endotoxins	HEK-Blue™ LPS Detection Kit 2	Assay for the detection and quantification of biologically active LPS	13
	HEK-Blue™ Selection	Antibiotics for maintenance of HEK-Blue™ cells	13
	HEK-Blue™ Endotoxin Standard	Standardized <i>E. coli</i> 055:B5 LPS	13
	QUANTI-Blue™ Solution	Alkaline phosphatase detection medium (liquid form)	13
Fungi	Fungin™	Reagent for preventing and eliminating fungal contamination	14

## REFERENCES

- Drexler H.G and Uphoff C.C., 2002. Mycoplasma contamination of cell cultures: incidence, sources, effects, detection, elimination, prevention. *Cytotechnology*. 39:75.
- Armstrong E.E. *et al.*, 2010. The scope of mycoplasma contamination within the biopharmaceutical industry. *Biologicals*. 38(2):211-3.
- Takeuchi O. *et al.*, 2000. Cutting Edge: Preferentially the R-Stereoisomer of the Mycoplasma Lipopeptide Macrophage-Activating Lipopeptide-2 Activates Immune Cells Through a Toll-Like Receptor 2- and MyD88- Dependent Signaling Pathway. *J. Immunol.* 164:554-557.
- Miller C.J. *et al.*, 2003. Mycoplasma infection significantly alters microarray gene expression profiles. *Biotechniques*. 35(4):812-4
- Rottem S. 2003. Interaction of Mycoplasmas with host cells. *Physiol Rev.* 83:417.
- Zakharova E. *et al.*, 2010. Mycoplasma suppression of THP-1 Cell TLR responses is corrected with antibiotics. *PLoS One*. 5(3):e9900
- Molla Kazemiha V *et al.*, 2011. Efficiency of Plasmocin™ on various mammalian cell lines infected by mollicutes in comparison with commonly used antibiotics in cell culture: a local experience. *Cytotechnology*. Dec; 63(6):609-20
- Uphoff CC *et al.*, 2012. Treatment of mycoplasma contamination in cell cultures with Plasmocin. *J Biomed Biotechnol.* 2012:26767.
- Rongvaux A *et al.*, 2014. Development and function of human innate immune cells in a humanized mouse model. *Nat Biotechnol.* 32(4):364-72.
- Baronti C. *et al.*, 2013. Mycoplasma removal: simple curative methods for viral supernatants. *J Virol Methods*. 187(2):234-7.
- Deng F. *et al.*, 2012. Generation of induced pluripotent stem cells from human Tenon's capsule fibroblasts. *Mol Vis.* 18:2871-81.
- Molla Kazemiha V. *et al.*, 2019. Effectiveness of Plasmocure in elimination of mycoplasma species from contaminated cell cultures: a comparative study versus other antibiotics. *Cell. J.* 21(2):143-149.
- Romorini L *et al.*, 2013. Effect of antibiotics against Mycoplasma sp. on human embryonic stem cells undifferentiated status, pluripotency, cell viability and growth. *Plos One*. Jul 30;8(7):e70267
- Ryan J., 2008. Understanding and managing cell culture contamination. Corning Life Sciences, technical literature.
- Lincoln C. & Gabridge M. 1998. Cell culture contamination: sources, consequences, prevention, and elimination. *Methods in cell biology*. 57:49-65.
- Fogh J., 1973. Contamination in Tissue Culture, published by Academic Press Inc. 4. McGarrity G.J. & Coriell LL., 1971. Procedures to reduce contamination of cell cultures. *In Vitro*. 6(4):257-65.
- McGarrity, G. J. *et al.*, 1984. Cytogenetic effects of mycoplasma infection of cell cultures: a review. *In vitro*. 20(1):1-18.
- Gray JS. *et al.*, 2010. Got black swimming dots in your cell culture? Identification of *Achromobacter* as a novel cell culture contaminant. *Biologicals*. 38(2):273-7
- McGowan JE Jr., 2006. Resistance in nonfermenting gram-negative bacteria: multidrug resistance to the maximum. *Am J Med.* 119(6 Suppl 1):S29-36; discussion S62-70.
- Gorbet MB. & Sefton MV., 2005. Endotoxin : the uninvited guest. *Biomaterials*. 26:6811-6817.
- Leider, R. *et al.*, 2013. Endotoxin - The invisible companion in biomaterial research. *Tissue Eng. Part B. Rev.* 19(5):391-402.
- Nims, RW. & Price P.J. 2017. Best practices for detecting and mitigating the risk of cell culture contaminants. *In vitro Cell. Dev. Biol. Anim.* 53(10):872-79.
- Rathinam VAK, *et al.*, 2019. Innate immunity to intracellular LPS. *Nat. Immunol.* 20(5): 527-533.
- Mather J. & Roberts E. 1998. Contamination: How to avoid it, recognize it, and get rid of it. In: *Introduction to cell and tissue culture: theory and technique*. Chapt 7. p99-9.

# Mycoplasma?

NEW



## MycoStrip™ Mycoplasma Detection Kit

*A new way to detect mycoplasma in your cell culture*

- ✓ **Simple** - No special equipment required
- ✓ **Rapid** - Performed in 1 hour, less than 15 min of hands-on time
- ✓ **Clear** - One band: negative, Two bands: positive for Mycoplasma
- ✓ **Specific** - No cross reactivity with eukaryotic or other bacterial DNA
- ✓ **Sensitive** - Able to detect as low as  $10^{-10^2}$  CFU/ml



**InvivoGen**  
www.invivogen.com

Europe Tel: +33 562 71 69 39  
USA Tel: +1 888 457 5873  
Asia Tel: +852 3622 3480

Fax: +33 562 71 69 30  
Fax: +1 858 457 5843  
Fax: +852 3622 3483

info.eu@invivogen.com  
info@invivogen.com  
info.hk@invivogen.com